

CONSUMER TESTIMONIALS

Congratulations on purchasing the FaceMaster®!

Included in this package you'll find everything you need to get amazing results. Before you dive into the Instructional DVD and Manual, I want to take you through the Quick Start so you can see and feel how this incredible product works immediately. Go ahead... take it out of the box right now!

QUICK START

1

Plug in wand until you hear a "click".

2

Fill well with FaceMaster Conductive Solution.

3

Place foam caps onto wand tips.

4

Dip wand tips into conductive solution.

5

Press power button.

6

Touch wand tips together to activate.

7

When activated, display screen will indicate a space over the eye area with a count down timer.

Ready to begin!

"I love the Facemaster®. I look like I had a facelift. I can't believe this product worked like you said it would. I have taken 10 years off my face in only a week. Thank you. Thank you. Thank you. I am 58 years old and had a bad sagging face. Not any more, thanks to you. I feel great.

-Doris Duran • age 58

"I can't believe how quickly this started to tone my face! As soon as my daughter saw the results, she got her own. The FaceMaster® is the best money I have ever spent on my face!"

- Barb Dreher • age 73

"The FaceMaster® is the easiest part of my beauty routine. In just a few minutes a day, my cheeks have firmed up, my eyelids look younger, and my friends are all asking for my secret!"

- Denise Snyder • age 42

If you have any questions at all, we have a friendly voice on the end of the telephone waiting for your call. We realize most of our customers have never tried anything like this. Feel free to call for any type of assistance or to order more supplies. We are here for you!



Made exclusively for: FaceMaster® of Beverly Hills • Calabasas • CA • 91302

1.800.770.2521 • facemaster.com

300021-A

suzanne
SOMERS

FaceMaster
facial toning system



QUICK START